


Fluorescence + Absorbance Spectrometer Fiber Optic



- Simultaneous obtained of all spectra
- Time Monitoring of spectra
- Integration time controlling using software
- Mathematical processing (Averaging, Sum, Smooth ...)
- Variable Excite Wavelength

طیف سنج فلورسانس + جذب

Model: FL+Abs-Arn-VF20



Fluorescence + Absorbance Spectrometer

FL+Abs-Arn-VF20

Specifications:

Optical Specifications

Wavelength Range	250-800 nm	Emission	FL
Wavelength Range	310-550 nm	Excite Variable	FL
Wavelength Range	200-800 nm		Abs
Resolution	< +1.0 nm		
Scan Time	2-1000 ms		
Typical Scan Time	0.1 s		
Stray Light	< 0.03 %		
Wavelength Accuracy	< ±0.1 nm		
Photometric Stability	< 0.01 A/min at 0 A , 500 nm		
Photometric Noise	< 0.01 A 80 scan at 0 A, 500 nm		
Wavelength Reproducibility	< ±0.02 nm for 10 scans (NIST 2034)		
Photometric Accuracy	< ± 0.01 A at 440.0 465.0 590.0 635.0 nm (NIST 930e)		
Baseline Flatness	< 0.01 A 0.5 second blank, 0.5 second scan at 0 A , 500 nm		
Light Source	Deuterium + Halogen/Tungsten Lamp or Xenon Lamp		
Detector	2000 Multichannel		



Physical Specifications

Communication	USB2
Windows	7 , 8 , 10
Software	Spectroscopy Station
Weight	6 Kg
HxWxD	22x22x28 cm
Power	220V AC 50/60 Hz

Applications:

